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
- 213189 Anti-Jk^a for card method Coombs-reactive, polyclonal, human 5 ml
 213184 Anti-Jk^b for card method Coombs-reactive, polyclonal, human 5 ml

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INTENDED USE

Polyclonal Coombs-reactive Anti-Jk^a and -Jk^b reagents are produced from human plasma that contains antibody of IgG-type, which reacts specific with the corresponding antigen. The reagents are used to in-vitro determine qualitative whether red blood cells possess or lack the corresponding blood group antigens Jk^a and Jk^b. The reagents are intended to be used by qualified and technical personnel only.

PRINCIPLE OF PROCEDURE

The test method used with these reagents is based on the principle of column agglutination. Normal human erythrocytes, possessing one of these antigens, will be recognized and coated by the corresponding specific antibody and then the cells will be agglutinated by a secondary antibody that reacts with human IgG-molecules.

REAGENTS

The listed reagents are available in following formulation:

Anti-Jk^a for card method Coombs-reactive, polyclonal, human

Anti-Jk^b for card method Coombs-reactive, polyclonal, human

The reagents contain <0.1% (w/v) sodium azide as preservative. Beside the parts active antibody and human serum, the reagents contain sodium chloride, macromolecules and bovine albumin, which has been tested and certified by the US Veterinary service inspectors.

WARNING

The reagents are prepared from human plasma. The raw material for these reagents are tested for HBsAg, HIV 1/2- and HCV-antibodies and found to be negative. Nevertheless as biological products it should be looked upon as potentially infectious because of never complete exclusion of danger through excitants of disease. The reagents contain sodium azide, that may be toxic and may react with lead or copper to form highly explosive salts.

On disposal, flush with large quantities of water.

For the reasons mentioned above the reagents should be handled with proper care.

STORAGE REQUIREMENT

Store opened and unopened products at 2 to 8°C.

May be at room temperature while in use.

In principle, store and use the reagents to declared expiry date only.

REMARKS

- With each testing positive and negative controls should be performed.
- Inappropriate storage impairs efficacy of the reagents.
- Weak turbidity of the reagents does not affect its reactivity. Bacteria and chemical contamination of the product should be avoided. If a visible change is detected, the reagents should no longer be used, this sign may indicate a microbiological contamination.
- Strength of positive reactions also depends on age of used blood.
- Use the Grifols card centrifuge. The use of another card-specific centrifuge (each card centrifuge has its specified unchangeable g-force) may lead to false results due to the changed g-force.
- The test method identified below is for manual testing only. When using automated or semi-automated instruments, follow the procedures that are contained in the operator's manual provided by the device manufacturer. Laboratories must follow approved validation procedures.
- For usage of these reagents all effective national laws, directives and guidelines have to be observed, in Germany especially the „Richtlinien zur Gewinnung von Blut und Blutbestandteilen und zur Anwendung von Blutprodukten (Hämotherapie)“¹ in its actual form.
- The information on the use of the Grifols test cards in the relevant insert must be observed.

SAMPLE PREPARATION

- Blood sample should be collected by approved medical procedure.
- Blood samples to be tested should be used as soon as possible after blood collection to reduce the risk of false-positive and false-negative results due to improper storage or contamination of the reagents.
If a delay in testing occurs, samples should be stored at 2 to 8°C.
Blood drawn into EDTA should be tested within 7 days and samples treated with sodium citrate within 14 days after collection.
Blood bag / Donor Blood can be tested until the expiry date.

REAGENT PREPARATION

There is no preparation of the reagents required. Take and use the reagents directly from the vials.

PROCEDURE

Not provided material additionally needed

Card Method

- Cards: Grifols „DG Gel Coombs“
- Microliter pipette
- Tubes
- Grifols Card centrifuge
- Grifols Card Incubator
- DG Gel Sol (card specific diluent)

Test procedure

Card Method

- Prepare 0,8 % suspension of red blood cells in DG Gel Sol (card specific diluent).
- Add 50 µL of appropriate cell suspension to each marked micro tube.
- Add 25 µL of appropriate reagent to each micro tube.
- Incubate card at 37°C for 15 minutes in a Grifols incubator.
- Centrifuge the card in a Grifols card centrifuge with the for this centrifuge unchangeable g-force.
- Check macroscopically for agglutination within 30 minutes.
- Document the results.

INTERPRETATION OF RESULTS

Positive results (+): visible agglutination of erythrocytes is a positive result and indicates the presence of the corresponding antigen.

Negative results (-): No visible agglutination of erythrocytes is a negative result and indicates the absence of the corresponding antigen.

Read and interpret the results of the card method according to the Grifols card instruction for use.

LIMITATIONS OF THE PROCEDURE

- Inaccuracy at compliance with instructions written under section “Procedures” and “Interpretation of results” may lead to incorrect results.
- No valid conclusion concerning the test result can be reached, if controls with uncertain or false results occur.
- Enzyme treated erythrocytes or addition of bovine albumin and/or other solutions containing protein may cause unspecific reactions.
- Hemolyzed, turbid, contaminated or clotted blood samples must not be used in this test.
- Due to variability of antigen expression on human red blood cells, reactivity of the reagents, mentioned above, against certain phenotypes may give weaker reactivity compared to control cells.
- No one specific antiserum or technique can be guaranteed to detect all rare, weak or variant antigens².
- Red blood cells coated with alloantibodies or autoantibodies of the same or similar specificity as that reagent used for the test (i.e. cells that are positive in the direct antiglobulin-test (DAT)) are not suitable for this test procedure.
- Red blood cells with a positive direct Coombs-test may cause false-positive reactions in card method.
- It is described in the literature, that samples from patients treated with anti-CD38 monoclonal antibodies can cause false positive results in the Coombs test³.
- Pay attention to all statements to limitations in inserts of Grifols card.

PERFORMANCE CHARACTERISTICS

A performance evaluation for the products was conducted.

Different samples material was used and compared with a reference method.

Anti-Jk^a agreement of the results and performance data:

When testing with the „DG Gel Coombs“ card Blood bag with the additive solution SAG-M and PAGGS-M, as well as EDTA- und Citratblood were used.

Method	Positive	Negative	Spensitivity	Specificity
Grifols Coombs Card	245/245	93/94	100%	98,94%

Anti-Jk^b agreement of the results and performance data:

When testing with the „DG Gel Coombs“ card Blood bag with the additive solution PAGGS-M, as well as EDTA- und Citratblood were used.

Method	Positive	Negative	Spensitivity	Specificity
Grifols Coombs Card	54/54	20/20	100%	100%

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PRESENTATION

213189 Anti-Jk^a for card method Coombs-reactive, polyclonal, human 5 ml
213184 Anti-Jk^b for card method Coombs-reactive, polyclonal, human 5 ml

730-15-5206 Version 006 / 18. May 2020

CE 0483

Antitoxin GmbH - Industriestrasse 88 - 69245 Bammental - Germany

Anti-Jk^a for card method Coombs-reactive, polyclonal, human

Anti-Jk^b for card method Coombs-reactive, polyclonal, human

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Antitoxin GmbH - Industriestraße 88 - 69245 Bammental - Nemecko