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Anti-CDE (RH2,1,3), monoclonal IgM + IgG (human) OrthoClone

REF 690141A

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REF 690141A

INTENDED USE

Monoclonal Anti-CDE test serum is produced from cell culture supernatants of heterohybridoma celllines. The cells are secreting an antibody of IgM- or IgG-type, which reacts specific with the corresponding blood group antigens. Each antibody is human protein. The test serum is used for In-Vitro-Diagnostic, to determine qualitative whether human red blood cells possess or lack the corresponding blood group antigen C, D and E. The test serum is intended to be used by qualified and technical personnel only.

PRINCIPLE OF PROCEDURE

The test methods used with this test serum are based on the principle of agglutination. Normal human erythrocytes, possessing one of the corresponding antigens, will agglutinate in the presence of the specific antibody directed towards the antigen.

REAGENT

The listed reagent contains antibodies of the following cell clones:

Anti-CDE	monoclonal, human IgM+IgG	clones: MS-24, MS-201, MS-80, MS-26
Anti-C	(clone: MS-24 (IgM))	
Anti-D	(clones: MS-201 (IgM), MS-26 (IgG))	
Anti-E	(clone: MS-80 (IgM))	

The test serum contains <0.1% (w/v) sodium azide as preservative. Additionally the test serum is comprised of active antibody, sodium chloride, macromolecules and bovine albumin, which has been tested and certified by the US Veterinary service inspectors.

CAUTION: This test serum is prepared from supernatants of cell cultures. As biological product it should be looked upon as potentially infectious because of never complete exclusion of danger through excitants of disease. The test serum contains sodium azide that may be toxic and may react with lead or copper to form highly explosive salts. On disposal, flush with large quantities of water. Found from the above reasons test serum should be handled with proper care.

STORAGE REQUIREMENTS

Store at +2 to +8°C (unopened/opened), or at room temperature while in use. Do not use reagent beyond their labeled expiry date.

REMARKS

1. It is recommended that each lot of test serum should be tested with appropriate positive and negative controls.
2. Inappropriate storage impairs efficacy of the test serum.
3. Weak turbidity of the reagent does not affect its reactivity. Bacteria and chemical contamination of the product should be avoided. If a visible change is detected, the reagent should no longer be used, this sign may indicate a microbiological contamination.
4. Strength of positive reactions also depends on age of used blood
5. Centrifugation outside the specified speed range may lead to false results.
6. The procedures identified below are for manual testing only. When using automated or semi-automated instruments, follow the procedures that are contained in the operator's manual provided by the device manufacturer. Laboratories must follow approved validation procedures.
7. For usage of this test serum all effective national laws, directives and guidelines must be followed. [in Germany especially the „Richtlinien zur Gewinnung von Blut und Blutbestandteilen und zur Anwendung von Blutprodukten (Hämotherapie)“¹] in its current version.

REAGENT PREPARATION

1. Blood sample should be collected by approved medical procedure.
2. Blood sample to be tested should be used as soon as possible after blood collection to reduce the risk of false-positive and false-negative results due to improper storage or contamination of the sample.
If a delay in testing occurs, sample should be stored at +2 to +8°C.
Blood drawn into EDTA should be tested within 7 days and sample treated with sodium citrate within 14 days after collection.
Blood bag / Donor Blood can be tested until the expiry date.

REAGENT PREPARATION

There is no specific preparation of the test serum required. Use test serum directly from the vial.

PROCEDURE

Material required but not provided:

Slide Method

1. Glass slide
2. Pasteur pipette
3. Mixing stick
4. Timer

Tube Centrifugation Method

1. Test tubes, 10 x 75 mm or 12 x 75 mm
2. Pipettes designed to deliver 100 µL, 50 µL
3. Timer
4. Centrifuge
5. Inkubator
6. Isotonic saline (0.85 – 0.9% sodium chloride)
7. Disposable pipette tips
8. Anti-Human-Globulin-Serum (Coombs-serum)

Test procedure

Slide Method

1. Use erythrocyte sediment only.
2. Place one drop (approximately 50 µL) of test serum on a marked glass slide.
3. Using a Pasteur pipette add one drop erythrocyte sediment (approximately 50 µL) to the glass slide.
4. Mix the erythrocytes with test serum well with a stick and spread to a circle with a diameter of approximately 2 cm.
5. By slightly rotating the slide, check for agglutination within 1 minute (reaction starts within seconds).
6. Document the result.

Tube Centrifugation Method

1. Prepare 2% to 5% suspensions of red blood cells in isotonic saline only (red blood cells may be washed 1-3 times with isotonic saline).
2. At first put 100 µL of the reagent in a marked tube subsequently add 100 µL of appropriate cell suspension in the tube.
Alternative one drop = approximately 50 µL cell suspension can be added to one drop = approximately 50 µL test serum.
3. Mix well by slightly shaking.
4. Incubate tube at room temperature for 15 min.
5. Centrifugation of tube for 1 minute at 2.000 rpm (approximately 800-1.000 g).
6. Gently resuspend the red cells and check macroscopically for agglutination within 3 minutes.
7. Document the result.
8. Incubate all tube with negative or weak result at +37 °C for 15 minutes.
9. Centrifugation of tube for 1 minute at 2.000 rpm (approximately 800-1.000 g).
10. Gently shake the red cells completely from the bottom of the tube and check macroscopically for agglutination within 3 minutes.
11. Document the result.
12. For all negative or weakly positive results, an indirect Coombs test must be performed.

Indirect Coombs-test

13. Repeat step 1.-3. with fresh material and incubate tube at +37 °C for 15 min or continue after reading of result directly.
14. Wash red cells 3 times with (cold) isotonic saline.
15. Add 100 µL Anti-Human-Globulin-Serum (Coombs-serum) to each tube.
16. Centrifugation of tube for 1 min at 1.000 rpm (approximately 180-270 g).
17. Gently shake the red cells completely from the bottom of the tube and check macroscopically for agglutination within 3 minutes.
18. Document the result.

INTERPRETATION OF RESULTS

"Slightly rotating / shaking" at Slide- and at Tube Centrifugation Method

Positive result (+): visible agglutination of erythrocytes indicates the presence of one or more of the corresponding antigens.

Negative result (-): no visible agglutination of erythrocytes indicates the absence of all of the corresponding antigens.

LIMITATIONS OF THE PROCEDURE

1. The "Procedures" and "Interpretation of Results" sections must be followed closely to assure the accuracy of the test result.
2. No valid conclusion concerning the test result can be reached, if controls with uncertain or false results occur.
3. Enzyme treated erythrocytes or addition of bovine albumin and/or other solutions containing protein may cause unspecific reactions.
4. Hemolyzed, turbid, contaminated or clotted samples should not be used.
5. Due to variability of antigen expression, reactivity of this test serum against certain phenotypes may give weaker reactivity compared to control cells.
6. Unspecific reactions might appear due to drying of the reaction-formation or if the slide is heated.
7. No one specific antiserum or technique can be guaranteed to detect all rare, weak or variant antigens ².
8. Red blood cells coated with alloantibodies or autoantibodies of the same or similar specificity as the test serum (i.e., cells that are positive in the direct antiglobulin-test (DAT)) are not suitable for this test procedure.
9. As described in the literature, samples from patients treated with anti-CD38 monoclonal antibodies can cause false positive results in the Coombs test.⁵
10. Newborn blood is excluded from testing.
11. The reagent anti-CDE does not react with rG cells.
12. The IgM Anti-D component will react negatively with D category VI cells, the IgG Anti-D component reacts positive with these cells by the indirect antiglobulin method.

PERFORMANCE CHARACTERISTICS

In compliance with Common Technical Specifications (CTS Commission Decision of 03. February 2009) a performance evaluation was conducted for the product. Different samples (donor, patient, panel blood) was used and compared with other products. The determined sensitivity and specificity are as follows

Product	Positive Blood	Negative Blood	False positive	False Negative	Sensitivity	Spezifitcity
Anti- CDE monoclonal, human IgM+IgG	1208	705	1	1	99,86%	99,92%

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LEGEND OF SYMBOLS

 Store from - to	 REF Product Code	 Expiration Date	 Manufacturer as to 98/79/EU
 LOT Lot	 CLON Clone(s)	 IVD In vitro diagnostic medical device	 CE EU CE-symbol
 UDI Unique Device Identification	 Consult Instrucion for use		

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Antitoxin GmbH, Industriestraße 88, 69245 Bammental, Deutschland