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 +49 (0) 6223/ 8661-0  +49 (0) 6223/ 8661-13  gara@antitoxin-gmbh.de

INSTRUCTIONS FOR USE (EN)

Anti-Kp^a, Coombsreactive, polyclonal (human)

REF 690652A

INTENDED USE

Poly clonal Coombsreactive Anti-Kp^a reagent is produced from human plasma that contain a specific antibody of IgG-type, which reacts exclusively with the corresponding antigen. The reagent is used for In-Vitro-Diagnostic, to determine whether red blood cells possess or lack the corresponding blood group antigens Kp^a. The reagent is intended to be used by qualified and technical personnel only.

PRINCIPLE OF PROCEDURE

The test method used with this reagent is based on the principle of agglutination. Normal human erythrocytes, possessing the appropriate antigen, will be recognized and coated with the corresponding specific antibody and then the cells will be agglutinated by a secondary antibody, which reacts with human IgG-molecules.

REAGENT

The listed reagent is available in following formulation:

Anti-Kp^a Coombs-reactive, polyclonal, human

The reagent contains <0.1% (w/v) sodium azide as preservative. Beside the parts, active antibody and the human serum, the reagent contain sodium chloride, macromolecules and bovine albumin, which has been tested and certified by the US Veterinary service inspectors.

CAUTION: This reagent is prepared from human plasma. The raw materials for this product have been tested for HBsAg, HIV 1/2- and HCV antibodies and found to be negative. Regardless, as biological product, it should be looked upon as potentially infectious because of never complete exclusion of danger through excipients of disease. The reagent contains sodium azide that may be toxic and may react with lead or copper to form highly explosive salts. On disposal, flush with large quantities of water.

For the reasons mentioned above, the reagent should be handled with proper care.

STORAGE REQUIREMENTS

Store at +2 to +8 °C (unopened / opened), or at room temperature while in use. Do not use reagent beyond its labelled expiration date!

REMARKS

1. It is recommended that each lot of reagent be tested with appropriate positive and negative controls.
2. Inappropriate storage impairs efficacy of the reagent.
3. Weak turbidity of the reagent does not affect its reactivity. Bacteria and chemical contamination of the product should be avoided. If a visible change is detected, the reagent should no longer be used, this sign may indicate a microbiological contamination.
4. Strength of positive reactions also depends on age of blood used.
5. Centrifugation outside the specified speed range may lead to false results.
6. The test method identified below is for manual testing only. When using automated or semi-automated instruments, follow the procedures that are contained in the operator's manual provided by the device manufacturer. Laboratories must follow approved validation procedures.
7. For usage of this reagent all effective national laws, directives and guidelines must be followed in its current version.
[In Germany especially: „Richtlinien zur Gewinnung von Blut und Blutbestandteilen und zur Anwendung von Blutprodukten (Hämotherapie)“¹.]

SAMPLE PREPARATION

1. Blood sample should be collected by approved medical procedure.
2. Blood sample to be tested should be used as soon as possible after blood collection to reduce the risk of false-positive and false-negative results due to improper storage or contamination of the sample.
If a delay in testing occurs, sample should be stored at +2 to +8°C.
Blood drawn into EDTA should be tested within 7 days and sample treated with sodium citrate within 14 days after collection.
Blood bag / Donor Blood can be tested until the expiry date.

REAGENT PREPARATION

There is no specific preparation of the reagent required. Use reagent directly from the vials.

PROCEDURE

Material required but not provided:

At Tube Centrifugation Method

1. Test tubes, 10 x 75 mm or 12 x 75 mm
2. Pipettes designed to deliver 50 µL/100 µL
3. Disposable pipette tips
4. Timer
5. Incubator
6. Centrifuge
7. Isotonic saline (0.85 – 0.9% sodium chloride)
8. Anti-Human-Globulin-serum (Coombs-serum)

Test procedure

Tube Centrifugation Method

1. Prepare a 2% to 5% suspension of red blood cells in isotonic saline (cells may be washed one time or up to three times with isotonic saline).
2. Add 100 µL (alternative: one drop, approximately 50 µL) of the reagent to a marked tube.
3. Add 100 µL (alternative: one drop, approximately 50 µL) of appropriate cell suspension to the tube.
4. Mix well by slightly shaking.
5. Incubate tubes at +37 °C for 30 min.
6. Wash red cells 3 times with (cold) isotonic saline.
7. Add 100 µL Anti-Human-Globulin-serum (Coombs-serum) to the tube, release the cells from the bottom of the tube by slightly shaking and mix with the Coombs-Serum / AHG-Serum.
8. Centrifugation of tube for 1 min at 1.000 rpm (approximately 180-270 g).
9. Gently shake the red cells from the bottom of the tube and check macroscopically for agglutination within 3 minutes.
10. Document the result.

INTERPRETATION OF RESULTS

"Slightly shaking" at Tube Centrifugation Method

Positive result (+): visible agglutination of erythrocytes indicates the presence of the corresponding antigen.

Negative result (-): no visible agglutination of erythrocytes indicates the absence of the corresponding antigen.


LIMITATIONS OF PROCEDURE

1. The "Procedure" and "Interpretation of Results" sections must be followed closely to assure the accuracy of the test results.
2. No valid conclusion concerning the test result can be obtained, if controls with unclear or false results should occur.
3. Enzyme treated erythrocytes or addition of bovine albumin and/or other solutions containing protein may cause unspecific reactions.
4. Hemolyzed, turbid, contaminated or clotted samples should not be used.
5. Due to variability of antigen expression, reactivity of these reagents against certain phenotypes may produce a weaker reaction compared to control cells.
6. No one specific antiserum or technique can be guaranteed to detect all rare, weak or variant Antigens.²
7. Red blood cells coated with alloantibodies or autoantibodies of the same or similar specificity as the reagent used for the test (i.e., cells that are positive in the direct antiglobulin-test [DAT]) are not suitable for this test procedure.
8. Red blood cells with a positive direct Coombs-test may cause false-positive reactions.
9. As described in the literature, samples from patients treated with anti-CD38 monoclonal antibodies can cause false positive results in the Coombs test.⁵

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SYMBOL - LEGENDE

| | | | |
|--|---|--|--|
|  Store from - to |  Product Code |  Expiration Date |  Manufacturer according to 98/79/EU |
|  Lot |  Clone(s) |  In vitro diagnostic medical device |  EU CE-symbol |
|  Unique Device Identification |  Consult Instruction for use | | |

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Antitoxin GmbH, Industriestraße 88, 69245 Bammental, Njemačka