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**Anti-Xg<sup>a</sup> Coombs-reactive, polyclonal, human**

For indirect Coombs-Test  
FOR IN VITRO DIAGNOSTIC USE ONLY

**INTENDED USE**

The reagent is used to in-vitro determine qualitative whether human red blood cells possess or lack the corresponding blood group antigen Xg<sup>a</sup>.

The reagent is intended to be used by qualified and technical personnel only to perform immunohematology screening tests as part of the practice of transfusion medicine in the general population.

The test method used with this reagent is based on the principle of agglutination, performed not automated.

Normal human erythrocytes, possessing the corresponding antigen, will be recognized and coated with the corresponding specific antibody and then the cells will be agglutinated by a secondary antibody, which reacts with human IgG-molecules..

**INDICATION / CONTRA-INDICATION**

The polyclonal Anti-Xg<sup>a</sup> blood grouping reagent is used to test patient or donor red cells for the presence of the Xg<sup>a</sup> antigen. Typing of donor cells facilitates the selection of suitable antigen-negative units for transfusion to patients with this antibody. Cell typing also serves as final verification of the identification of Anti-Xg<sup>a</sup> in patient or donor sera.

There is no contra-indication to perform the in-vitro-test on blood samples.

The product was validated with sample collected in Europe from patients of unknown ethnic background.

The approximate frequencies of Xg<sup>a</sup> antigen:

Phenotype	Women	Men
positive	89%	66%
negative	11%	34%

**REAGENTS**

The listed reagents are available in following formulation:

Anti-Xg<sup>a</sup> Coombs-reactive polyclonal, human

The reagent is produced from human plasma that contains a specific antibody of IgG-type, which reacts exclusively with the corresponding antigen.

The antibodies are human protein.

The reagent contains <0.1% (w/v) sodium azide as preservative.

In addition to the active antibody component the test serum contains, human serum, sodium chloride, macromolecules, bovine albumin, which has been tested negative for Vesicular Stomatitis Virus and Bluetongue.

The BSA was derived from US sourced animals from USDA and APHIS approved facilities, to be used for in vitro diagnostic reagents, in compliance with EC 1069/2009 and EU 142/2011.

**WARNING**

This reagent is prepared from human plasma. Nevertheless, The raw materials for these products have been tested for HBsAg, HIV 1/2- and HCV antibodies and found to be negative. As biological product, it should be looked upon as potentially infectious because of never complete exclusion of danger through excipients of disease. The reagent contains sodium azide that may be toxic and may react with lead or copper to form highly explosive salts.

On disposal, flush with large quantities of water.

For the reasons mentioned above the reagent should be handled with proper care.

**STORAGE REQUIREMENT**

Store opened and unopened products at +2 to +8°C. May be at room temperature while in use.

For the stimulation of an In Use Study, the sera were stored 30 times at room temperature for 2 hours and subsequently showed no differences in the qualitative tests until the expiry date.

In principle, store and use the reagent to declared expiry date only.

(Format for the expiry date: Year xxxx Month xx Day xx)

**REMARKS**

1. With each testing positive and negative controls should be performed.
2. Test the Anti-Human Globulin Serum (Coombs serum / AHG serum) for functionality with positive and negative controls before use.
3. Inappropriate storage impairs efficacy of the reagent.
4. Weak turbidity of the reagent does not affect its reactivity. Bacteria and chemical contamination of the product should be avoided. If a visible change (increase in turbidity or colour change due to temperature) is detected, the reagent should no longer be used, this sign may indicate a microbiological contamination.
5. Do not use leaking, unlabeled or broken vials.
6. Strength of positive reactions also depends on age of used blood.
7. Centrifugation outside the specified speed range may lead to false results.
8. The test method identified below is for manual testing only and must be carried out according to the instructions for use.
  - a) In case of changes in technology / deviations from the Instructions for Use
  - b) use of automated or semi-automated systems
 the laboratories have to follow the procedures that are contained in the operator's manual provided by the device manufacturer and carry out validations according to recognized procedures.
9. For usage of the reagent all effective national laws, directives and guidelines have to be observed in its actual form, in Germany especially the „Richtlinien zur Gewinnung von Blut und Blutbestandteilen und zur Anwendung von Blutprodukten (Hämotherapie)“<sup>1</sup>

**SAMPLE PREPARATION**

1. Blood sample should be obtained using an acceptable phlebotomy technique. Sample drawn into tubes containing EDTA, Sodium Citrate, CPD-A, ACD or Blood bag containing PAGGS-M may be used.
2. Blood samples to be tested should be used immediately after blood collection to reduce the risk of false-positive and false-negative results due to improper storage or contamination of the samples. Samples that cannot be tested immediately should be stored at +2 to +8°C. Blood drawn into EDTA should be tested within 7 days and samples treated with sodium citrate, CPD-A and ACD within 14 days after collection. Blood bag / Donor Blood (with PAGGS-M) can be tested until the expiry date.
3. Do not freeze samples.

**REAGENT PREPARATION**

There is no preparation of the reagent required. Take and use the reagent directly from the vials.

**ADDITIONAL REQUIRED MATERIAL NOT SUPPLIED AND THE ASSOCIATED PROCEDURAL INSTRUCTION****Tube-Centrifugation-Method:****Materials:**

1. Tubes (10 x 75 mm or 12 x 75 mm)
2. Microliter pipette
3. Centrifuge
4. Incubator
5. Timer
6. isotonic saline (0,85 - 0,9% sodium chlorid)
7. Anti-Human-Globulin-serum (Coombs-serum)

**Procedure:**

1. Prepare 2 to 5 % suspensions of red blood cells in isotonic saline (red blood cells may be washed 1 - 3 times with isotonic saline).
2. At first put 100 µL of appropriate reagent in a marked tube subsequently add 100 µL of appropriate cell suspension in the tube. Alternative one drop = approximately 50 µL cell suspension can be added one drop = approximately 50 µL test serum.
3. Mix Erythrocytes-/Reagent mixture well by slightly shaking.
4. Incubate tubes for 30 minutes at +37 °C in an incubator.
5. Wash red cells 3 times with (cold) isotonic saline.
6. After the third wash, aspirate the NaCl supernatant.
7. Then add 100 µL Anti-Human-Globulin-serum (Coombs-serum) to the tube, gently shake o remove the cell button from the bottom of the tube and mix with the Coombs serum / AHG serum.
8. Centrifugation of tube for 1 minute at 1.000 rpm (approximately 180-270 x g).
9. Gently shake the red cells completely from the bottom of the tube and check macroscopically for agglutination within 3 minutes.
10. Document the result.

**INTERPRETATION OF RESULTS**

Positive results (+): If agglutination of the erythrocytes occurs within the accepted limitations of the test procedure, the test result is positive and indicates the presence of the corresponding antigen.

Negative results (-): If agglutination of the erythrocytes does not occur within the accepted limitations of the test procedure, the test result is negative and indicates the absence of the corresponding antigen.

The reading and interpretation of the results after "careful shaking" at the Tube Centrifugation Method:

Negative	No detectable agglutinates, homogeneous red coloration of the liquid.
Positive	One complete agglutinate.
	No complete agglutinate, some individual agglutinates
	Red colouration of the liquid containing only small / miniature agglutinates.

**LIMITATIONS OF THE PROCEDURE**

1. Inaccuracy at compliance with instructions written under section "Procedures" and "Interpretation of results" may lead to incorrect results.
2. No valid conclusion concerning the test result can be reached, if controls with uncertain or false results occur.
3. Enzyme treated erythrocytes or addition of bovine albumin and / or other solutions containing protein may cause unspecific reactions.
4. Hemolyzed, turbid, contaminated or clotted blood samples must not be used in this test.
5. A red blood cell suspension with a concentration that deviates from the indicated concentration may lead to false positive or false negative results.
6. The addition of volumes that deviate from the volumes specified in the method may lead to altered reaction behavior.
7. Due to variability of antigen expression on human red blood cells, reactivity of the reagent mentioned above, against certain phenotypes, may give weaker reactivity compared to control cells
8. Red blood cells coated with alloantibodies or autoantibodies of the same or similar specificity as the appropriate reagent (i.e., cells that are positive in the direct antiglobulin-test [DAT]) are not suitable for this test procedure.
9. Red blood cells with a strong positive direct Coombs-test may give false positive results in rare cases.
10. No one specific antiserum or technique can be guaranteed to detect all rare, weak or variants of antigens.<sup>2</sup>
11. It is described in the literature, that samples from patients treated with anti-CD38 monoclonal antibodies can cause false positive results in the Coombs test.<sup>3</sup>

**INCIDENTS RELATED TO THE DEVICE**

Any serious incident that has occurred in relation to the device must be reported to the manufacturer and the competent authority of the Member State where the user and/or patient is established.

**PERFORMANCE CHARACTERISTICS**

A performance evaluation for the products was carried out.  
The required samples were used and compared with other reference methods / products.

Technique	Anti-Xg <sup>®</sup> Coombs-reactive, polyclonal, human			
	Positive Bleeds n	Sensitivity	Negative Bleeds n	Specificity
Tube methode	15 / 15	100 %	110 / 110	100 %

Diagnostic Sensitivity: The probability that the device gives a positive result in the presence of the target marker.

Diagnostic Specificity: The probability that the device gives a negative result in the absence of the target marker.

Anti-Xg<sup>®</sup> Coombs-reactive, polyclonal, human is equivalent and does not differ in quality from comparable reagents available on the market.

**DIFFERENCES BETWEEN BATCHES**

Validation between three batches over the entire shelf life showed no differences.

**INTERFERENCE STUDY**

The interference studies showed no impairment for the qualitative test when using the following interfering substances:

Heparin 720 U/dl, Albumin 15000 mg/dl, Triglycerides 1500 mg/dl, Bilirubin 40mg/dl, Ethanol 620 mg/dl, Glucose 1000 mg/dl.


For the anticoagulants (EDTA, Sodium Citrate, ACD, CPD-A, PAGGS-M) three times the recommended concentration was tested.

**SUMMARY OF SAFETY AND PERFORMANCE**

The Summary of Safety and Performance of this reagent is available via ANTITOXIN ([www.antitoxin-gmbh.de](http://www.antitoxin-gmbh.de)) and can be accessed via the EUDAMED database.

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 Product code	 Batch
 Store from - to	 Expiration Date
 In-Vitro Diagnostic	 EU CE symbol
 Manufacture accoding to (EU) 2017/746	 Consult instruction for use
 Unique Device Identification	 Distributor

**REF**213081 Anti-Xg<sup>®</sup> Coombs-reactive, polyclonal, human 2 ml**CE 0123**
 ANTITOXIN GmbH Industriestraße 88 69245 Bammental Germany

 +49 (0) 6223/ 8661-0

 +49 (0) 6223/ 8661-13

 [gara@antitoxin-gmbh.de](mailto:gara@antitoxin-gmbh.de)
 02.063- / Version R001 / 2024-04-10

Marking of changes

Underlined: addition or substantial change; ◆ Deletion of texts

# Anti-Xg<sup>a</sup> Coombs-reactive polyclonal, human

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ANTITOXIN GmbH Industriestraße 88 69245 Bammental Alemanha