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In-vitro-Diagnostikum

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REF 690221A

INTENDED USE

Monoclonal agglutinating Anti-C^w reagent is produced from cell culture supernatants of a hetrohybridoma-cell line. The cells secrete antibodies of IgM-type, which reacts specific with the corresponding blood group antigen. The antibody is human protein. The reagent is used to in-vitro determine qualitative whether human red blood cells possess or lack the blood group antigen C^w. The reagent is intended to be used by qualified and technical personnel only.

PRINCIPLE OF PROCEDURE

The test methods used with this reagent is based on the principle of agglutination. Normal human erythrocytes, processing the corresponding antigen, will agglutinate in the presence of the specific antibody.

REAGENT

The listed reagent contains antibodies of the following clone:

Anti-C^w monoclonal, human IgM clone: MS-110

The reagent contains <0.1% (w/v) sodium azide as preservative. Beside the part active antibody, the reagent contains sodium chloride, macromolecules and bovine albumin, which has been tested and certified by the US Veterinary service inspectors.

CAUTION: The reagent is prepared from supernatants of cell cultures. Regardless, as biological product, it should be looked upon as potentially infectious because of never complete exclusion of danger through excitants of disease. The reagent contains sodium azide that may be toxic and may react with lead or copper to form highly explosive salts. On disposal, flush with large quantities of water. For the reasons mentioned above, the reagent should be handled with proper care.

STORAGE REQUIREMENTS

Store at +2 to +8 °C (unopened / opened) or at room temperature while in use. Do not use reagent beyond its labelled expiration date.

REMARKS

1. With each testing positive and negative controls should be performed.
2. Inappropriate storage impairs efficacy of the reagent.
3. Weak turbidity of the reagent does not affect its reactivity. Bacteria and chemical contamination of the product should be avoided. If a visible change is detected, the reagent should no longer be used, this sign may indicate a microbiological contamination.
4. Strength of positive reactions also depends on age of used blood
5. Centrifugation outside the specified speed range may lead to false results.
6. The test methods identified below are for manual testing only. When using automated or semi-automated instruments, follow the procedures that are contained in the operator's manual provided by the device manufacturer. Laboratories must follow approved validation procedures.
7. For usage of the reagent all effective national laws, directives and guidelines in its current version have to be observed, in Germany especially the "Richtlinien zur Gewinnung von Blut und Blutbestandteilen und zur Anwendung von Blutprodukten (Hämotherapie)"¹.

SAMPLE PREPARATION

1. Blood sample should be collected by approved medical procedure.
2. Blood sample to be tested should be used as soon as possible after blood collection to reduce the risk of false-positive and false-negative results due to improper storage or contamination of the sample. If a delay in testing occurs, sample should be stored at +2 to +8°C. Blood drawn into EDTA should be tested within 7 days and sample treated with sodium citrate within 14 days after collection. Blood bag / Donor Blood can be tested until the expiry date.

REAGENT PREPARATION

There is no specific preparation of the reagent required. Take and use reagent directly from the vial.

PROCEDURE

Not provided material, additionally needed:

Slide Method	Tube Centrifugation Method
1. glass slide	1. tubes, 10 x 75 mm or 12 x 75 mm
2. pasteur pipette	2. pipettes designed to deliver 50 µL/100 µL
3. mixing stick	3. disposable pipette tips
4. timer	4. timer
	5. centrifuge
	6. isotonic saline (0.85 – 0.9% sodium chloride)

Test procedure

Slide Method

1. Use erythrocyte sediment only.
2. Place one drop (approximately 50 µL) of the reagent on a marked glass slide.
3. Using a Pasteur pipette add one drop erythrocyte sediment (approximately 50 µL) to the drop of reagent on the glass slide.
4. Mix the erythrocytes with reagent well with a stick and spread to a circle with a diameter of approximately 2 cm.
5. By slightly rotating the slide, check for agglutination within 1 minute (reaction starts within seconds).
6. Document the result.

Tube Centrifugation Method

1. Prepare 2% to 5% suspensions of red blood cells in isotonic saline (red blood cells may be washed 1-3 times with isotonic saline).
2. Add 100 µL (alternative: one drop, approximately 50 µL) of reagent to a marked tube.
3. Add 100 µL (alternative: one drop, approximately 50 µL) of appropriate cell suspension to the tube.
4. Mix well by shaking slightly.
5. Incubate tube at room temperature for 15 min.
6. Centrifuge tube for 1 minute at 2.000 rpm (approximately 800–1.000 x g).
7. Gently shake the red cells from the bottom of the tube and check macroscopically for agglutination within 3 minutes.
8. Document the result.

INTERPRETATION OF RESULTS

"Slightly rotation" at Slide Method / "slightly shaking" at Tube Centrifugation Method:

Positive result (+): visible agglutination of erythrocytes indicates the presence of the corresponding antigen.

Negative result (-): no visible agglutination of erythrocytes indicates the absence of the corresponding antigen.


LIMITATION OF THE PROCEDURE

1. The "Procedure" and "Interpretation of Results" sections must be followed closely to assure the accuracy of the test results.
2. No valid conclusion concerning the test result can be obtained, if controls with unclear or false results should occur.
3. Enzyme treated erythrocytes or addition of bovine albumin and/or other solutions containing protein may cause unspecific reactions.
4. Hemolyzed, turbid, contaminated or clotted samples should not be used.
5. With the slide method, unspecific reactions might appear due to drying of the reaction-formation or if the slide is heated.
6. Due to variability of antigen expression on human red blood cells, reactivity of the reagent against certain phenotypes, may give weaker reactivity compared to control cells.
7. No one specific antiserum or technique can be guaranteed to detect all rare, weak or variant Antigens.²
8. Red blood cells coated with alloantibodies or autoantibodies of the same or similar specificity as that reagent used for the test (i.e., cells that are positive in the direct antiglobulin-test (DAT)) are not suitable for this test procedure.

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SYMBOL - LEGENDE

 Store from - to	REF Product Code	 Expiration Date	 Manufacturer according to 98/79/EU
LOT Lot	CLON Clone(s)	IVD In vitro diagnostic medical device	 EU CE-symbol
UDI Unique Device Identification	 Consult Instrucion for use		

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Antitoxin GmbH, Industriestraße 88, 69245 Bammental, Deutschland