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qara@antitoxin-gmbh.de

For Slide, Tube, Card and Microplate Method
FOR IN VITRO DIAGNOSTIC USE ONLY

INTENDED USE

Monoclonal Anti-C, -c, -E, -e- reagents are produced from cell culture supernatants of heterohybridoma cell lines. The cells are secreting an antibody of IgM-type that reacts specific with the corresponding blood group antigen. The antibody in each case is human protein. The reagents are used for In-Vitro-Diagnostic, to determine qualitative, whether human red blood cells possess or lack the corresponding blood group antigens. C, c, E, and e.

The testsera are intended to be used by qualified and technical personnel only.

PRINCIPLE OF PROCEDURE

The test method used with these reagents are based on the principle of agglutination. Normal human erythrocytes, possessing the corresponding antigen, will agglutinate in the presence of the specific antibody directed toward the antigen.

REAGENTS

The listed monoclonal, human IgM reagents are available in two different variations. The differences depend on cell clones only:

- Anti-C clones: MS-24, P3x25513G8 Anti-C clone: MS-273
- Anti-c clone: MS-33 Anti-c clone: MS-35
- Ank-E clones: MS-258, 906 Anti-E clones: MS-12, MS-260
- Anti-e clones: MS-16, MS-21, MS-63 Anti-e clones: MS-62, MS-69

The reagents contain <0.1% (w/v) sodium azide as preservative. Additionally the reagents are prepared of active antibody, sodium chloride, macromolecules and bovine albumin, which has been tested and certified by the US Veterinary service inspectors.

WARNING

These reagents are prepared from supernatants of cell cultures. As biological products it should be looked upon as potentially infectious because of never complete exclusion of danger through excitants of disease. The reagents contain sodium azide that may be toxic and may react with lead or copper to form highly explosive salts.

On disposal, flush with large quantities of water.

For the reasons mentioned above, reagent should be handled with proper care.

STORAGE REQUIREMENT

Store opened and unopened products at 2 to 8°C. May be at room temperature while in use. In principle, store and use the reagents to declared expiry date only.

REMARKS

1. With each testing positive and negative controls should be performed.
2. Inappropriate storage impairs efficacy of the reagents.
3. Weak turbidity of the reagents mentioned above does not affect its reactivity. Bacteria and chemical contamination of the products should be avoided. If a visible change is detected, the reagents should no longer be used, this sign may indicate a microbiological contamination.
4. Strength of positive reactions also depends on age of used blood
5. Centrifugation outside the specified speed range may lead to false results. With the card method, the use of another card-specific centrifuge (each card centrifuge has its specified unchangeable g-force) may lead to false results due to the changed g-force
6. The test methods identified below are for manual testing only. When using automated or semi-automated instruments, follow the procedures that are contained in the operator's manual provided by the device manufacturer. Laboratories must follow approved validation procedures.
7. For usage of this testsera all effective national laws, directives and guidelines have to be observed, in Germany especially the „Richtlinien zur Gewinnung von Blut und Blutbestandteilen und zur Anwendung von Blutprodukten (Hämotherapie)“ * in its current version.
8. The information on the use of the test cards in the relevant insert must be observed.

SAMPLE PREPARATION

1. Blood sample should be collected by approved medical procedure.
2. Blood sample to be tested should be used as soon as possible after blood collection to reduce the risk of false-positive and false-negative results due to improper storage or contamination of the sample. If a delay in testing occurs, samples should be stored at 2 to 8°C. Blood drawn into EDTA should be tested within 7 days and samples treated with sodium citrate within 14 days after collection. Blood bag / Donor Blood can be tested until the expiry date.

REAGENT PREPARATION

There is no preparation of the reagents required. Use reagents directly from the vials.

PROCEDURE

Not provided material, additionally needed

- Slide Method: glass slide; Pasteur pipette; mixing stick; timer tubes (10 x 75 mm or 12 x 75 mm); microliter pipette; timer; centrifuge; incubator; isotonic saline (0,85 - 0,9% sodium chloride)
- Tube Centrifugation Method: Grifols „DG Gel Neutral“; tubes; microliter pipette; centrifuge; timer; corresponding card centrifuge; card specific diluent; isotonic saline (0,85 - 0,9% sodium chloride)
- Card Method: microplate with 96 U-wells; microliter pipette; timer; centrifuge; microplate shaker; microplate centrifuge; isotonic saline (0,85 - 0,9% sodium chloride)
- Microplate Method: microplate with 96 U-wells; microliter pipette; timer; centrifuge; microplate shaker; microplate centrifuge; isotonic saline (0,85 - 0,9% sodium chloride)

Test procedure

Slide Method

1. Use erythrocyte sediment only.
2. Place one drop (approximately 50 µL) of appropriate reagent on a marked glass slide.
3. Using a Pasteur pipette add one drop erythrocyte sediment (approximately 50 µL) to the drop of reagent on the glass slide.
4. Mix the erythrocytes with reagent well with a stick and spread to a circle with a diameter of approximately 2 cm.
5. By slightly rotating the slide, check for agglutination within 1 minute (reaction starts within seconds).
6. Document the result.

Tube Centrifugation Method

1. Prepare a 2% to 5% suspension of red blood cells in isotonic saline (red blood cells may be washed 1-3 times with isotonic saline).
2. At first put 100 µL of the appropriate reagent in a marked tube subsequently add 100 µL of appropriate cell suspension in the tube. Alternative one drop = approximately 50 µL cell suspension can be added to one drop = approximately 50 µL test serum.
3. Mix Erythrocytes-/Reagent mixture well by slightly shaking.
4. Incubate tube at room temperature for 15 min.
5. Centrifugation of tube for 1 minute at 2.000 rpm (approximately 800 - 1.000 x g).
6. Gently shake the red cells completely from the bottom of the tube and check macroscopically for agglutination within 3 minutes.
7. Document the result.

Card Method

1. Prepare 0,8 % suspension of red blood cells in card specific diluent (red blood cells may be washed 1-3 times with isotonic saline).
2. Add 50 µL of appropriate cell suspension to a marked micro tube.
3. Add 25 µL of appropriate testserum to the micro tube.
4. Centrifuge the card in appropriate card centrifuge with the for this centrifuge unchangeable g-force after latest 30 minutes.
5. Check macroscopically for agglutination within 30 minutes.
6. Document the result.

Microplate Method

1. Prepare 2-5% suspension of red blood cells in isotonic saline (red blood cells may be washed 1-3 times with isotonic saline).
2. Add 50 µL of appropriate reagent to a marked well.
3. Add 50 µL of appropriate cell suspension to the well.
4. Shake microplate for 30 seconds on shaker with medium speed.
5. Centrifugation of the microplate in appropriate microplate centrifuge for 30 seconds at 400 x g.
6. Shake microplate for 30 seconds on a microplate shaker with medium speed.
7. Evaluated macroscopically for agglutination directly after centrifugation.
8. Document the result.
9. Tests with negative or doubtful results have to be incubated for 5 to 10 minutes at room temperature.
10. Repeat steps 5 to 8 after incubation at room temperature.

INTERPRETATION OF RESULTS

"Slightly rotating/shaking"

at Slide -/ Tube Centrifugation - and Microplate Method:

Positive results (+): visible agglutination of erythrocytes is a positive result and indicates the presence of the corresponding antigen.

Negative results (-): No visible agglutination of erythrocytes is a negative result and indicates the absence of the corresponding antigen.

Read and interpret the results of the gel card according to the instruction for use.

LIMITATIONS OF THE PROCEDURE

1. Inaccuracy at compliance with instructions written under section "Procedures" and "Interpretation of results" may lead to incorrect results.
2. No valid conclusion concerning the test result can be reached, if controls with uncertain or false results occur.
3. Enzyme treated erythrocytes or addition of bovine albumin and/or other solutions containing
4. Hemolyzed, turbid, contaminated or clotted samples should not be used.
5. With the slide method, unspecific reactions might appear due to drying of the reaction formation or if the slide is heated.
6. Due to variability of antigen expression, reactivity of these reagents against certain phenotypes may give weaker reactivity compared to control cells.
7. Some red blood cells may have developed RH variants that may not be recognized with monoclonal antibodies of a clone. Therefore the tests should / have to be carried out with at least two different clones / cloning mixtures as required in some countries. Deviating test results require further testing.
8. No one specific antiserum or technique can be guaranteed to detect all rare, weak or variant antigens.²
9. Red blood cells coated with alloantibodies or autoantibodies of the same or similar specificity as the reagent (i.e., cells that are positive in the direct antiglobulin-test (DAT)) are not suitable for this test procedure.
10. Red blood cells coated with antibodies (cells that are positive in the direct antiglobulin-test (DAT)) may give false-positive results in card method. These cells react positive without testserum too.
11. Pay attention to all statements to limitations in inserts of used gel cards.
12. Negative test results with anti-e (clones: MS-62, MS-69) in microplate method have to be confirmed by an other technique, as false negative results may occur.
13. Newborn blood is excluded for the microtiter plate test method
14. With cord cells anti-c (clone: MS-33) may possess anti-I/i cold agglutinin activity and therefore show false positive results.

PERFORMANCE CHARACTERISTICS

In compliance with Common Technical Specifications (CTS Commission Decision of 03. February 2009) a performance evaluation was conducted. Different samples (donor,patient, panel blood) was used and compared with other products

Product	Positive	Negative	False Positive	False Negative
Anti-C (MS-24, P3x25513G8)	649	305	0	0
Anti-C (MS-273)	649	305	0	0
Anti-c (MS-33)	777	177	0	0
Anti-c (MS-35)	777	177	0	0
Anti-E (MS-3258, 906)	249	705	0	0
Anti-E (MS-12, MS-260)	249	705	0	0
Anti-e (MS-16, MS-21, MS-63)	938	16	0	0
Anti-e (MS-62, MS-69)	1138	63	0	0*

The calculated values found for:

Product	Sensitivity	Specificity
Anti-C (MS-24, P3x25513G8)	100%	100%
Anti-C (MS-273)	100%	100%
Anti-c (MS-33)	100%	100%
Anti-c (MS-35)	100%	100%
Anti-E (MS-3258, 906)	100%	100%
Anti-E (MS-12, MS-260)	100%	100%
Anti-e (MS-16, MS-21, MS-63)	100%	100%
Anti-e (MS-62, MS-69)	100%*	100%

* for Anti-e (clones MS-62, MS-69) in microplate method 42 samples were tested false negative, so in this method sensitivity was 97,12% only.

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ANTITOXIN GmbH Industriestraße 88 69245 Bammental Alemanha