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### Anti-K (KEL1), monoclonal IgM (human) OrthoClone

**REF** 690461A

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**REF** 690461A

#### INTENDED USE

Monoclonal agglutinating Anti-K (KEL1) reagent is produced from cell culture supernatants of a hetero-hybridoma-cell line. The cells secrete an antibody of IgM-type which reacts specific with the corresponding blood group antigen. The antibody is human protein. The reagent is used for In-Vitro-Diagnostic, to determine qualitative whether human red blood cells possess or lack the corresponding blood group antigen K. The procedures described are for the listed manual methods and application with Ortho BioVue® Reverse Diluent cassettes on the ORTHO VISION™ Analyzer and ORTHO AutoVue® Innova.

The test serum is intended to be used by qualified technical personnel only.

#### PRINCIPLE OF PROCEDURE

The procedures used with this reagent is based on the principle of agglutination. Normal human erythrocytes, possessing the corresponding antigen, will agglutinate in the presence of the specific antibody directed towards the antigen.

#### REAGENT

The listed reagent contain antibody of the following clone:

Anti-K(KEL1) monoclonal, human IgM clone: AEK-4

This reagent contains <0.1% (w/v) sodium azide as preservative. Additionally the reagent is prepared of active antibody, sodium chloride, macromolecules and bovine albumin, which has been tested and certified by the US Veterinary service inspectors.

**CAUTION:** This reagent is prepared from supernatants of cell cultures. As biological product it should be looked upon as potentially infectious because of never complete exclusion of danger through excitants of disease. The reagent contain sodium azide, which may be toxic and may react with lead or copper to form highly explosive salts.

On disposal, flush with large quantities of water. For the reasons mentioned above please handle reagent with proper care.

#### STORAGE REQUIREMENTS

Store at +2 to +8 °C (unopened / opened), or shortly at room temperature while in use. If the devices ORTHO VISION™ Analyzer or ORTHO Auto Vue® Innova are used, the reagent may remain a maximum of 5 days on the device (+15 bis +30 °C). Do not use reagent beyond the stated expiration date!

#### REMARKS

1. It is recommended that each lot of reagent be tested with appropriate positive and negative controls.
2. Inappropriate storage impairs efficacy of the reagent.
3. Weak turbidity of the reagent does not affect its reactivity. Bacteria and chemical contamination of the product should be avoided. If a visible change is detected, this sign may indicate a microbiological contamination, the reagent should no longer be use.
4. Strength of positive reactions also depends on age of blood used.
5. Centrifugation outside the specified speed range may lead to false results. With the card method, the use of another card-specific centrifuge (each card centrifuge has its specified unchangeable g-force) may lead to false results due to the changed g-force.
6. Blood samples to be tested should be used as soon as possible. If a delay in testing occurs, samples should be stored at +2 to +8 °C. Blood drawn into EDTA should be tested within 7 days and sample drawn into natrium citrate within 14 days after collection. Blood bag / Donor Blood can be tested until the expiry date.
7. The procedure identified below is for manual testing only. When using automated or semi-automated instruments, follow the procedures that are contained in the operator's manual provided by the device manufacturer. Laboratories must follow approved validation procedures.
8. For usage of this reagent all effective national laws, directives and guidelines must be followed. In Germany especially: „Richtlinien zur Gewinnung von Blut und Blutbestandteilen und zur Anwendung von Blutprodukten (Hämotherapie)“<sup>1</sup> in its current version.
9. The information on the use of the Ortho BioVue® Reverse Diluent Test-Kassette in the associated instructions must be observed.

#### REAGENT PREPARATION

There is no specific preparation of the reagent required. Take and use reagent directly from the vial.

#### PROCEDURE

Material required but not provided:

Tube Centrifugation Method

1. Test tubes, 10 x 75 mm or 12 x 75 mm
2. Pipettes designed to deliver 50µL/100 µL
3. Timer
4. Centrifuge
5. Isotonic saline (0,85 - 0,9% sodium chloride)
6. Disposable pipette tips

Microplate Method

1. Microplate with U-bottom, (optional: pre-treated)
2. Centrifuge, suitable for centrifugation of microplates
3. Microplate-carrier for centrifuge (option)
4. Microplate shaker
5. Microplatetest mirror (option)
6. Test tubes, 10 x 75 mm or 12 x 75 mm
7. Pipette for 50 µL
8. Timer
9. Isotonic saline (0.85 - 0.9 % sodium chloride)
10. Disposable pipette tips

Application of BioVue® Cassettes

1. Ortho BioVue® Reverse Diluent Cassettes
2. Isotonic saline (0,85 - 0,9 % sodium chloride)
3. Test tubes, 10 x 75 mm or 12 x 75 mm
4. Mikroliter pipette for 10 µL, 40 µL and 50 µL
5. Disposable pipette tips
6. Timer
7. Centrifuge
8. Ortho BioVue® System Work Rack
9. Ortho BioVue® System centrifuge

#### Test procedure

##### Tube Centrifugation Method

1. Prepare 2% to 5% suspensions of red blood cells in isotonic saline only (red blood cells may be washed one time or up to three times with isotonic saline).
2. At first put 100 µL of appropriate reagent in a marked tube subsequently, add 100 µL of appropriate cell suspension in the tube. Alternative one drop = approximately 50 µL cell suspension can be added to one drop = approximately 50 µL test serum.
3. Mix well by slightly shaking.
4. Incubate tube at room temperature for 15 min.
5. Centrifugation of tube for 1 minute at 2.000 rpm (approximately 800-1.000 g).
6. Gently resuspend the red cells and within 3 minutes check macroscopically for agglutination.
7. Document the result.

##### Microplate Method

1. Prepare 2% to 5% suspensions of red blood cells in isotonic saline (red blood cells may be washed one time or up to three times with isotonic saline).
2. Add one drop (approximately 50 µL) of appropriate reagent into a marked well.
3. Add one drop (approximately 50 µL) of cell suspension into the well.
4. Mix microplate on a microplate shaker.
 

**Note:** Recommended:

  - (1) To mix: 10 to 30 seconds at medium strength.
  - (2) To resuspend: 10 to 30 seconds at medium strength respectively with duration and strength, allowing a full resuspension of sediment without destruction of positive reactions.
5. Centrifugation of microplate.
 

Recommended duration of centrifugation: 30 seconds at approximately 400 x g, or as long as necessary to produce best positive reactions with antigen positive cells, with used centrifuge and microplate, and still have an easy resuspension of antigen negative cells.

**Note:** The used centrifugal force should be the minimum force to produce a sediment and a clear supernatant. Because of multitude offered centrifuges a General speed or duration for best results cannot be recommended. Each laboratory has to test and validate its equipment for speed and duration to obtain best results.
6. Resuspend the erythrocytes buttons with the shaker and position microplate for reading.
7. Read microplate and document the result.
8. Tests with negative or doubtful results have to be incubated for 5 to 10 minutes at room temperature.
10. Repeat steps 5 to 7 after incubation at room temperature.
 

**Note:** All tests should be resuspended and evaluated directly after centrifugation.

##### Manual Test Method in Ortho BioVue® Reverse Diluent Cassettes

The procedure described below is for manual Ortho BioVue® Reverse Diluent Cassettes testing only. Please refer to the instructions for use of the Ortho BioVue® Reverse Diluent Cassettes, sections STORAGE REQUIREMENT, PRECAUTIONS, SPECIMEN COLLECTION and PREPARATION, REAGENT PREPARATION and INTERPRETATION OF RESULTS, as well as -where applicable- LIMITATIONS OF THE PROCEDURE.

1. For preparation of the cassettes, please refer to the instructions for use for Ortho BioVue® Reverse Diluent Cassettes, chapter „PROCEDURE“ point 1 and 2.
2. Prepare 3 % to 5% suspensions of red blood cells in isotonic sodium chloride solution. (red blood cells may be washed one time or up to three times with isotonic saline)
3. Add 10 µL of a appropriate 3 to 5% red blood cell suspension to the marked reaction chamber of the cassette.
4. Add 40 µL appropriate reagent to the appropriate reaction chambers of the cassette.  
(The pipette tips shall not contact the rim of the reaction chambers. If this should adversely happen, please replace the tip with a new one, before you will fill the next reaction chamber.)

5. Centrifuge the cassette in the Ortho BioVue® System centrifuge with the for this centrifuge unchangeable g-force.  
**Note:** Centrifugation must be performed within 30 minutes after addition of specimens to the reaction chambers.
6. After completion of test, read each column from the cassette's front and rear sides for agglutination and/or hemolysis.
7. Protocol the final test result from the cassette side with the more clearly expressed positive reading.

**Test in Ortho BioVue® Reverse Diluent Cassettes on ORTHO VISION™ Analyzer or ORTHO AutoVue® Innova**

Please follow the instructions in the user manual of the respective analyzer.

### INTERPRETATION OF RESULTS

"Slightly shaking" at Tube Centrifugation Method and using a microplate shaker at Microplate Method

**Positive result (+):** visible agglutination of erythrocytes indicates the presence of the corresponding antigen.

**Negative result (-):** no visible agglutination of erythrocytes indicates the absence of the corresponding antigen.

**Cassettes:** As for the manual method, for interpretation of results, please refer to the instructions for use of Ortho BioVue® Reverse Diluent Cassettes, chapter „INTER-PRATATION OF RESULTS“.  
 As for analyzers, reading and interpretation is effected automatically by each device.

### LIMITATIONS OF THE PROCEDURE

1. The "Procedure" and "Interpretation of Results" sections must be followed closely to assure the accuracy of the test results.
2. No valid conclusion concerning the test result can be reached, if controls with uncertain or false results occur.
3. Enzyme treated erythrocytes or addition of bovine albumin and/or other solutions containing protein may cause unspecific reactions.
4. Hemolyzed, turbid, contaminated or clotted blood samples must not be used in this test.
5. Due to variability of antigen expression, reactivity of this reagent against certain phenotypes may produce a weaker reaction compared to control cells.
6. No one specific antiserum or technique can be guaranteed to detect all rare, weak or variant antigens <sup>2</sup>.
7. Red blood cells coated with alloantibodies or autoantibodies of the same or similar specificity as the reagent (i.e., cells that are positive in the direct antiglobulin-test [DAT]) may not suitable for this test procedure.
8. Red blood cells coated with antibodies (cells that are positive in the direct antiglobulin-test (DAT)) may give false-positive results in cassettes method. These cells react positive without reagent too.
9. Pay attention to all statements to limitations in inserts of used Ortho BioVue® Reverse Diluent Cassettes and the user manual of respective analyzer.

### PERFORMANCE

In compliance with Common Technical Specifications (CTS Commission Decision of 03. February 2009) a performance evaluation was conducted. Different samples (donor, patient, neonates-, panel blood) was used and compared with other reference methods / products.

	Number of specimens tested	Specimens pos with Anti-K	Specificity	Specimens neg with Anti-K	Sensitivity
Tube / Mikroplate	310	196	100%	114	100%
Manual BioVue®	142	42	100%	100	99,47%
AutoVue®	102	16	100%	86	100%

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### LEGEND OF SYMBOLS

 Store from - to	 REF	Product Code	 Expiration Date	 Manufacturer as to 98/79/EU
 LOT	 i	Observe instructions for use	 IVD	 CE EU CE-symbol

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Antitoxin GmbH, Industriestraße 88, 69245 Bammental, Niemcy